



Sterile
Container System







Sterile Container System

GUARANTY

All products carrying the registered trademark and logo of BAHADIR are guaranteed for two years in respect of material and workmanship. We will replace or repair any damaged instrument free of charge during this period.

All rights reserved.

No part of this catalog may be reproduced in any form without written permission from BAHADIR

Suitability of our sterile containers for steam sterilization procedures¹⁾ as per DIN 58953, part 9 (issue May 1987)

Sterilization container	Packing		Steam sterilization procedure acc. to DIN 58946, part I				
	Description	Symbolic representation and packing system	Gravitation procedure	Pre - vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injection procedure
with filters (StB - F)	Textiles in sterilization container with filter in lid and bottom acc. to DIN 58952 part I		+	+	+	+	+
			0	0	+	+	+
	Textiles in sterilization container with filter in lid, bottom not perforated acc. to DIN 58952 part I		-	-	+	+	+
			-	-	+	+	0
	Textiles in sterilization container with filter in bottom lid, not perforated acc. to DIN 58952 part I		0	+	+	+	+
			-	0	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III in sterilization container with filter in lid and bottom acc. to DIN 58952 part I		+ ³⁾	+	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid and bottom not perforated acc. to DIN 58952 part I		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in bottom acc. to DIN 58952 part I		0	+	+	+	+
	with valves (StB - V)	Textiles in sterilization container with valves in lid and bottom acc. to DIN 58952 part I		-	0	+	+
			-	-	+	+	0
Textiles in sterilization container with valves in lid no valves on the bottom			-	-	+	+	-
			-	-	+	+	-
Instruments in wire basket in sterilization container with valve in lid and bottom			-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
Instruments in wire basket in sterilization container with valve in lid and not perforated			-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾

Explanation of symbols : + applicable - not applicable 0 limited application

1) In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures.

2) Only sterilization containers guaranteeing evaporation or discharge of condensate are suitable above.

3) Gravitation procedures having no drying device do not guarantee a sufficient after - drying effect.

Benefits, Sterilizing Procedures	6 - 7 - 8 - 9 - 10 - 11
1/1 Size Containers and Accessories	13 - 14 - 15 - 16 - 17
3/4 Size Containers and Accessories	19 - 20 - 21 - 22 - 23
1/2 Size Containers and Accessories	25 - 26 - 27 - 28 - 29
Wide Body Containers	31 - 32 - 33
Other Containers and Accessories	35 - 45
Flat Containers and Accessories	36 - 37
Mini Containers and Accessories	38 - 39
Dental Containers and Accessories	40
Implant Containers	41
Scopy Containers	41
Removal Containers and Accessories	42 - 43
Brush Containers	44 - 45
General Accessories	47 - 48 - 49
Trolleys	50
Numerical Index	51 - 52 - 53 - 54

Features & Benefits

- Long lasting sterilization.
- Maximum flexibility in the preparation of surgical sets and sterile material.
- Protects surgical instruments from damage.
- Standard size allow easy stacking of containers.
- Provides safe transportation of surgical instruments.
- Risk free return to storage against contamination.
- Cost-effective.

All dimensions comply with DIN 58952 standard,thus guaranteeing an economical rate of utilization of the sterilizers.

Container lids and bottoms are made of anodized aluminum-alloy,therefore,due to the thermal conductivity of the material an optimum drying phase is achived.

All standart size containers can be stacked,thus space and economic storage problems are eliminated.

Wire baskets can also be stacked to prevent valuable instruments from being damaged during preparation of the sets.

Removal containers; ideal for safe return of dirty instruments and textiles without any risk of contamination during transportation after operation.

Wrapping drapes provide a safe sterile medium in the containers and wire baskets; and guarantee the sterilized contents against contamination when they are taken out for operation.

Ensures an economic and long lasting sterilization.



Sterile Container System



We recommend following steam sterilization procedures for our containers with perforated lids and non-perforated bottoms, as specified in DIN 58946/6 standard.

For fractional flow procedure:

A few jets of saturated steam of $p_{\text{e}} \geq 1$ bar is generally sufficient to remove the air out of the sterilizing chamber.

For fractional vacuum procedure:

Repeated evacuation down to a pressure of ≈ 130 mbar in exchange with steam inflow up to a pressure above or below an atmospheric pressure.

Maintaining of steam inflow until the working overpressure is reached.

For steam injection procedure (vacuum procedure)

Single evacuation down to a pressure of ≈ 70 mbar with simultaneous low steam inflow.

Maintaining of steam inflow until the working pressure is reached.

For further assistance regarding all other gravitation and prevacuum procedures as specified in the DIN 58953/9 standard, please refer to the table on page 4.

basic model



Basic model containers are produced with perforated or non-perforated bases. Filters are not protected against contamination from liquids, dust and other mechanical effects. Lids are available in 6 different colours.





safe model

Safe containers are fitted with a safety lid protecting filters. They are supplied either with filters or without filters. Safety lids are supplied in 6 different colours.





General front view



Each container contains two identification labels; one for the surgical department which the container belongs to; one for the name of the operation. They are easily slid into their built-in slots and firmly secured by a spring enforcement.



The indicator labels are filled in by an authorized surgeon before the sterilization process and inserted into their spring-enforced slots.



The security seal is used to keep the container lock safe from unauthorized opening, and automatically breaks up only when the locking flap is lifted.




1-The locking flaps are movable 180° on a pivot and can be secured in any position.
 2-Silicon seals are provided around the lid and in the centre and periphery of the round filter retainer, completely recessed.

The filter retainer is locked by pressing down with a thumb in the direction of arrow "1" and open when pressed with a thumb in the direction of arrow "2".



To lift off the container lid; push the lock in the direction of arrow "1", and holding the lid latch, pull in the direction of arrow "2".
 For mounting the safety lid, simply place it on the container lid and press down.





1/1 size containers



Bottom Non Perforated



-  Y 110.10
-  Y 120.10
-  Y 130.10
-  Y 140.10
-  Y 150.10
-  Y 160.10

weight 2634 gr



outside dim.	580 x 280 x 100 mm
inside dim.	550 x 265 x 84 mm

Bottom Perforated



-  Y 111.10
-  Y 121.10
-  Y 131.10
-  Y 141.10
-  Y 151.10
-  Y 161.10

weight 3156 gr

-  Y 110.13
-  Y 120.13
-  Y 130.13
-  Y 140.13
-  Y 150.13
-  Y 160.13

weight 3129,5 gr



outside dim.	580 x 280 x 135 mm
inside dim.	550 x 265 x 115 mm

-  Y 111.13
-  Y 121.13
-  Y 131.13
-  Y 141.13
-  Y 151.13
-  Y 161.13

weight 3244 gr

-  Y 110.15
-  Y 120.15
-  Y 130.15
-  Y 140.15
-  Y 150.15
-  Y 160.15

weight 3219 gr



outside dim.	580 x 280 x 150 mm
inside dim.	550 x 265 x 130 mm

-  Y 111.15
-  Y 121.15
-  Y 131.15
-  Y 141.15
-  Y 151.15
-  Y 161.15

weight 3382 gr

-  Y 110.20
-  Y 120.20
-  Y 130.20
-  Y 140.20
-  Y 150.20
-  Y 160.20

weight 3427,5 gr



outside dim.	580 x 280 x 200 mm
inside dim.	550 x 265 x 180 mm

-  Y 111.20
-  Y 121.20
-  Y 131.20
-  Y 141.20
-  Y 151.20
-  Y 161.20

weight 3738,5 gr

-  Y 110.26
-  Y 120.26
-  Y 130.26
-  Y 140.26
-  Y 150.26
-  Y 160.26

weight 3919 gr



outside dim.	580 x 280 x 260 mm
inside dim.	550 x 265 x 240 mm

-  Y 111.26
-  Y 121.26
-  Y 131.26
-  Y 141.26
-  Y 151.26
-  Y 161.26

weight 4272,5 gr

safe model
Types and Dimensions

1/1 size containers

Bottom Non Perforated



-  Y 112.10
-  Y 122.10
-  Y 132.10
-  Y 142.10
-  Y 152.10
-  Y 162.10

weight 3486 gr



outside dim.	580 x 280 x 100 mm
inside dim.	550 x 265 x 84 mm

Bottom Perforated



-  Y 113.10
-  Y 123.10
-  Y 133.10
-  Y 143.10
-  Y 153.10
-  Y 163.10

weight 4007.5 gr

-  Y 112.13
-  Y 122.13
-  Y 132.13
-  Y 142.13
-  Y 152.13
-  Y 162.13







weight 3981.5 gr



outside dim.	580 x 280 x 135 mm
inside dim.	550 x 265 x 115 mm

-  Y 113.13
-  Y 123.13
-  Y 133.13
-  Y 143.13
-  Y 153.13
-  Y 163.13







weight 4095.5 gr

-  Y 112.15
-  Y 122.15
-  Y 132.15
-  Y 142.15
-  Y 152.15
-  Y 162.15

weight 4070.5 gr



outside dim.	580 x 280 x 150 mm
inside dim.	550 x 265 x 130 mm

-  Y 113.15
-  Y 123.15
-  Y 133.15
-  Y 143.15
-  Y 153.15
-  Y 163.15

weight 4233 gr

-  Y 112.20
-  Y 122.20
-  Y 132.20
-  Y 142.20
-  Y 152.20
-  Y 162.20


weight 4280 gr



outside dim.	580 x 280 x 200 mm
inside dim.	550 x 265 x 180 mm

-  Y 113.20
-  Y 123.20
-  Y 133.20
-  Y 143.20
-  Y 153.20
-  Y 163.20






weight 4589 gr

-  Y 112.26
-  Y 122.26
-  Y 132.26
-  Y 142.26
-  Y 152.26
-  Y 162.26

weight 4770.5 gr



outside dim.	580 x 280 x 260 mm
inside dim.	550 x 265 x 240 mm

-  Y 113.26
-  Y 123.26
-  Y 133.26
-  Y 143.26
-  Y 153.26
-  Y 163.26

weight 5140 gr

bottoms and lids

Bottom Non Perforated



Bottom Perforated



	outside dim.	inside dim.	weight
Y 105.10	580 x 280 x 100 mm	550 x 262 x 84 mm	1523 gr
Y 105.13	580 x 280 x 135 mm	550 x 262 x 117 mm	2019 gr
Y 105.15	580 x 280 x 150 mm	550 x 262 x 132 mm	2108 gr
Y 105.20	580 x 280 x 200 mm	550 x 262 x 182 mm	2317 gr
Y 105.26	580 x 280 x 260 mm	550 x 262 x 244 mm	2808 gr

	outside dim.	inside dim.	weight
Y 106.10	580 x 280 x 100 mm	550 x 262 x 84 mm	2044 gr
Y 106.13	580 x 280 x 135 mm	550 x 262 x 117 mm	2132 gr
Y 106.15	580 x 280 x 150 mm	550 x 262 x 132 mm	2270 gr
Y 106.20	580 x 280 x 200 mm	550 x 262 x 182 mm	2627 gr
Y 106.26	580 x 280 x 260 mm	550 x 262 x 244 mm	2851 gr

Types and Dimensions

580 x 280 mm

Lid Perforated

Y 100.01

Y 100.02

Y 100.03

Y 100.04

Y 100.05

Y 100.06

weight | 1111 gr



Safety Lid

Y 103.01

Y 103.02

Y 103.03

Y 103.04

Y 103.05

Y 103.06

weight | 851.5 gr



Special Accessories
Wire Baskets


	outside dim.	weight
Y 190.03	485 x 255 x 30 mm	1108 gr
Y 190.05	485 x 255 x 50 mm	1112 gr
Y 190.07	485 x 255 x 70 mm	1176 gr
Y 190.10	485 x 255 x 100 mm	1289 gr




	outside dim.	weight
Y 191.03	540 x 255 x 30 mm	1149 gr
Y 191.05	540 x 255 x 50 mm	1189 gr
Y 191.07	540 x 255 x 70 mm	1259 gr
Y 191.10	540 x 255 x 100 mm	1354 gr

Special Accessories
Wire Base



Y 193.00	545 x 255 x 30 mm	729 gr
----------	-------------------	--------

Drape Retainer


Y 195.00	545 x 255 mm
----------	--------------

Drape


Y 095.07	130 x 60 cm
Y 095.08	140 x 100 cm
Y 095.09	185 x 45 cm
Y 095.10	185 x 60 cm



3/4 size containers



Types and Dimensions

Bottom Non Perforated



-  Y 210.10
-  Y 220.10
-  Y 230.10
-  Y 240.10
-  Y 250.10
-  Y 260.10

weight 2268,5 gr



outside dim.	465 x 280 x 100 mm
inside dim.	425 x 260 x 84 mm

Bottom Perforated



-  Y 211.10
-  Y 221.10
-  Y 231.10
-  Y 241.10
-  Y 251.10
-  Y 261.10

weight 2503 gr

-  Y 210.13
-  Y 220.13
-  Y 230.13
-  Y 240.13
-  Y 250.13
-  Y 260.13

weight 2549,5 gr



outside dim.	465 x 280 x 135 mm
inside dim.	425 x 260 x 115 mm

-  Y 211.13
-  Y 221.13
-  Y 231.13
-  Y 241.13
-  Y 251.13
-  Y 261.13

weight 2715 gr

-  Y 210.15
-  Y 220.15
-  Y 230.15
-  Y 240.15
-  Y 250.15
-  Y 260.15

weight 2620,5 gr



outside dim.	465 x 280 x 150 mm
inside dim.	425 x 260 x 130 mm

-  Y 211.15
-  Y 221.15
-  Y 231.15
-  Y 241.15
-  Y 251.15
-  Y 261.15

weight 2718 gr

Types and Dimensions

Bottom Non Perforated



-  Y 212.10
-  Y 222.10
-  Y 232.10
-  Y 242.10
-  Y 252.10
-  Y 262.10

weight 2980 gr



outside dim.	465 x 280 x 100 mm
inside dim.	425 x 260 x 84 mm

Bottom Perforated



-  Y 213.10
-  Y 223.10
-  Y 233.10
-  Y 243.10
-  Y 253.10
-  Y 263.10

weight 3213.5 gr

-  Y 212.13
-  Y 222.13
-  Y 232.13
-  Y 242.13
-  Y 252.13
-  Y 262.13

weight 3260 gr



outside dim.	465 x 280 x 135 mm
inside dim.	425 x 260 x 115 mm

-  Y 213.13
-  Y 223.13
-  Y 233.13
-  Y 243.13
-  Y 253.13
-  Y 263.13


weight 3425 gr

-  Y 212.15
-  Y 222.15
-  Y 232.15
-  Y 242.15
-  Y 252.15
-  Y 262.15

weight 3975 gr




outside dim.	465 x 280 x 150 mm
inside dim.	425 x 260 x 130 mm

-  Y 213.15
-  Y 223.15
-  Y 233.15
-  Y 243.15
-  Y 253.15
-  Y 263.15

weight 3561 gr

Bottom Non Perforated



bottoms and lids

Bottom Perforated




	outside dim.	inside dim.	weight
Y 205.10	465 x 280 x 100 mm	425 x 262 x 85 mm	1396 gr
Y 205.13	465 x 280 x 135 mm	425 x 262 x 121 mm	1676 gr
Y 205.15	465 x 280 x 150 mm	425 x 262 x 135 mm	1792 gr

	outside dim.	inside dim.	weight
Y 206.10	465 x 280 x 100 mm	425 x 262 x 85 mm	1629 gr
Y 206.13	465 x 280 x 135 mm	425 x 262 x 121 mm	1842 gr
Y 206.15	465 x 280 x 150 mm	425 x 262 x 135 mm	1978 gr

Types and Dimensions

Lid Perforated




Safety Lid



Wire Baskets



	outside dim.	weight
Y 290.03	405 x 255 x 30 mm	972.5 gr
Y 290.05	405 x 255 x 50 mm	988 gr
Y 290.07	405 x 255 x 70 mm	1061 gr
Y 290.10	405 x 255 x 100 mm	1143 gr

Wire Base



Y 293.03	405 x 255 x 30 mm	598 gr
----------	-------------------	--------

Drape Retainer


405 x 255 mm Y 295.00

Drape


140 x 100 cm Y 095.08








1/2 size containers



Types and Dimensions

Bottom Non Perforated



-  Y 310.10
-  Y 320.10
-  Y 330.10
-  Y 340.10
-  Y 350.10
-  Y 360.10

weight | 1744,5 gr



outside dim.	285 x 280 x 100 mm
inside dim.	262 x 260 x 82 mm

Bottom Perforated



-  Y 311.10
-  Y 321.10
-  Y 331.10
-  Y 341.10
-  Y 351.10
-  Y 361.10

weight | 1890,5 gr

-  Y 310.13
-  Y 320.13
-  Y 330.13
-  Y 340.13
-  Y 350.13
-  Y 360.13

weight | 1947 gr



outside dim.	285 x 280 x 135 mm
inside dim.	262 x 260 x 116 mm

-  Y 311.13
-  Y 321.13
-  Y 331.13
-  Y 341.13
-  Y 351.13
-  Y 361.13






weight | 2125 gr

-  Y 310.15
-  Y 320.15
-  Y 330.15
-  Y 340.15
-  Y 350.15
-  Y 360.15

weight | 1981 gr



outside dim.	285 x 280 x 150 mm
inside dim.	262 x 260 x 132 mm

-  Y 311.15
-  Y 321.15
-  Y 331.15
-  Y 341.15
-  Y 351.15
-  Y 361.15

weight | 2462,5 gr

-  Y 310.20
-  Y 320.20
-  Y 330.20
-  Y 340.20
-  Y 350.20
-  Y 360.20







weight | 2101 gr



outside dim.	285 x 280 x 200 mm
inside dim.	262 x 260 x 182 mm

-  Y 311.20
-  Y 321.20
-  Y 331.20
-  Y 341.20
-  Y 351.20
-  Y 361.20







weight | 2513,5 gr

-  Y 310.26
-  Y 320.26
-  Y 330.26
-  Y 340.26
-  Y 350.26
-  Y 360.26

weight | 2483,5 gr



outside dim.	285 x 280 x 260 mm
inside dim.	262 x 260 x 238 mm


-  Y 311.26
-  Y 321.26
-  Y 331.26
-  Y 341.26
-  Y 351.26
-  Y 361.26

weight | 2657,5 gr

Types and Dimensions

Bottom Non Perforated









-  Y 312.10
-  Y 322.10
-  Y 332.10
-  Y 342.10
-  Y 352.10
-  Y 362.10

weight | 2162,5 gr


outside dim.	285 x 280 x 100 mm
inside dim.	262 x 260 x 82 mm

Bottom Perforated



-  Y 313.10
-  Y 323.10
-  Y 333.10
-  Y 343.10
-  Y 353.10
-  Y 363.10

weight | 2309,5 gr



-  Y 312.13
-  Y 322.13
-  Y 332.13
-  Y 342.13
-  Y 352.13
-  Y 362.13

weight | 2351 gr


outside dim.	285 x 280 x 135 mm
inside dim.	262 x 260 x 116 mm







-  Y 313.13
-  Y 323.13
-  Y 333.13
-  Y 343.13
-  Y 353.13
-  Y 363.13

weight | 2466 gr

-  Y 312.15
-  Y 322.15
-  Y 332.15
-  Y 342.15
-  Y 352.15
-  Y 362.15

weight | 2394 gr


outside dim.	285 x 280 x 150 mm
inside dim.	262 x 260 x 132 mm

-  Y 313.15
-  Y 323.15
-  Y 333.15
-  Y 343.15
-  Y 353.15
-  Y 363.15

weight | 2583 gr







-  Y 312.20
-  Y 322.20
-  Y 332.20
-  Y 342.20
-  Y 352.20
-  Y 362.20

weight | 2560 gr


outside dim.	285 x 280 x 200 mm
inside dim.	262 x 260 x 182 mm







-  Y 313.20
-  Y 323.20
-  Y 333.20
-  Y 343.20
-  Y 353.20
-  Y 363.20

weight | 2785,5 gr

-  Y 312.26
-  Y 322.26
-  Y 332.26
-  Y 342.26
-  Y 352.26
-  Y 362.26

weight | 2930 gr


outside dim.	285 x 280 x 260 mm
inside dim.	262 x 260 x 238 mm

-  Y 313.26
-  Y 323.26
-  Y 333.26
-  Y 343.26
-  Y 353.26
-  Y 363.26

weight | 3069,5 gr

bottoms and lids

Bottom Non Perforated



Bottom Perforated



	outside dim.	inside dim.	weight
Y 305.10	285 x 280 x 100 mm	262 x 262 x 82 mm	1096 gr
Y 305.13	285 x 280 x 135 mm	262 x 262 x 120 mm	1284 gr
Y 305.15	285 x 280 x 150 mm	262 x 262 x 134 mm	1327 gr
Y 305.20	285x 280 x 200 mm	262x 262 x 185 mm	1447 gr
Y 305.26	285 x 280 x 260 mm	262 x 262 x 240 mm	1817 gr

	outside dim.	inside dim.	weight
Y 306.10	285 x 280 x 100 mm	262 x 262 x 82 mm	1242 gr
Y 306.13	285 x 280 x 135 mm	262 x 262 x 120 mm	1398 gr
Y 306.15	285 x 280 x 150 mm	262 x 262 x 134 mm	1471 gr
Y 306.20	285x 280 x 200 mm	262x 262 x 185 mm	1672 gr
Y 306.26	285 x 280 x 260 mm	262 x 262 x 240 mm	2000 gr

Types and Dimensions

Lid Perforated



Y 300.01

Y 300.02

Y 300.03

Y 300.04

Y 300.05

Y 300.06

weight | 654 gr

285x 280 mm

Safety Lid



Y 303.01

Y 303.02

Y 303.03

Y 303.04

Y 303.05

Y 303.06

weight | 458 gr

285 x 280 mm


Special Accessories

Wire Baskets



	outside dim.	weight
Y 390.03	255 x 245 x 30 mm	637,5 gr
Y 390.05	255 x 245 x 50 mm	648,5 gr
Y 390.07	255 x 245 x 70 mm	691 gr
Y 390.10	255 x 245 x 100 mm	758,5 gr

Wire Base



	outside dim.	weight
Y 393.00	258 x 258 x 30 mm	351 gr

Special Accessories

Drape Retainer




255 x 255 mm Y 395.00

Drape



90 x 60 cm	Y 095.04
90 x 80 cm	Y 095.05
130 x 45 cm	Y 095.06
130 x 60 cm	Y 095.07



wide
body

containers





Wide body containers

Basic Model Types and Dimensions

Bottom Non Perforated



-  Y 110.62
-  Y 120.62
-  Y 130.62
-  Y 140.62
-  Y 150.62
-  Y 160.62







weight 3,857 gr



outside dim.	600 x 400 x 120 mm
inside dim.	550 x 350 x 100 mm

Bottom Perforated



-  Y 111.62
-  Y 121.62
-  Y 131.62
-  Y 141.62
-  Y 151.62
-  Y 161.62

weight 4101 gr



outside dim.	600 x 400 x 180 mm
inside dim.	550 x 350 x 160 mm

-  Y 110.68
-  Y 120.68
-  Y 130.68
-  Y 140.68
-  Y 150.68
-  Y 160.68

weight 4043,5 gr

-  Y 111.68
-  Y 121.68
-  Y 131.68
-  Y 141.68
-  Y 151.68
-  Y 161.68

weight 4581,5 gr

Special Accessories Wire Baskets



Y 190.55

470 X 337 x 51 mm



Y 191.60

533 X 337 X 101 mm

BAHADIR® Wide body containers

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



	outside dim.	inside dim.	weight
Y 105.62	600 X 400 X 120 mm	550 x 350 x 100 mm	2.189 gr
Y 105.68	600 X 400 X 180 mm	550 x 350 x 160 mm	2375 gr

	outside dim.	inside dim.	weight
Y 106.62	600 X 400 X 120 mm	550 x 350 x 100 mm	2.431 gr
Y 106.68	600 X 400 X 180 mm	550 x 350 x 160 mm	2911 gr

Types and Dimensions

600 X 400 mm

Lid Perforated



Y 100.51

Y 100.52


Y 100.53

Y 100.54


Y 100.55

Y 100.56

weight 1.666 gr



other containers









Flat containers

Types and Dimensions

Bottom Non Perforated



-  Y 310.05
-  Y 320.05
-  Y 330.05
-  Y 340.05
-  Y 350.05
-  Y 360.05

weight 1265 gr

Bottom Perforated



outside dim.	285 x 280 x 55 mm
inside dim.	270 x 270 x 33 mm



outside dim.	285 x 280 x 85 mm
inside dim.	270 x 270 x 64 mm

-  Y 310.08
-  Y 320.08
-  Y 330.08
-  Y 340.08
-  Y 350.08
-  Y 360.08


weight 1377 gr

-  Y 311.08
-  Y 321.08
-  Y 331.08
-  Y 341.08
-  Y 351.08
-  Y 361.08

weight 1529 gr




Bottom Non Perforated



Bottom Perforated



Lid Perforated



	outside dim.	inside dim.	weight
Y 305.05	285 x 280 x 55 mm	270 x 270 x 33 mm	611 gr
Y 305.08	285 x 280 x 85 mm	270 x 270 x 64 mm	703 gr

Y 306.08	285 x 280 x 85 mm	270 x 270 x 64 mm	857 gr
----------	-------------------	-------------------	--------


Special Accessories

Wire Baskets



255 x 245 x 30 mm Y 390.03
255 x 245 x 50 mm Y 390.05

Wire Base



Y 393.00	245 x 255 x 30 mm	351 gr
----------	-------------------	--------

Silicone Mat



245 x 245 mm Y 851.00

Drape Retainer



255 x 255 mm Y 395.00

Drape









90 x 60 cm Y 095.04
90 x 80 cm Y 095.05

Mini containers

Types and Dimensions

Bottom Non Perforated




-  Y 410.04
-  Y 420.04
-  Y 430.04
-  Y 440.04
-  Y 450.04
-  Y 460.04

weight 510.5 gr

Bottom Perforated



-  Y 411.04
-  Y 421.04
-  Y 431.04
-  Y 441.04
-  Y 451.04
-  Y 461.04

weight 589 gr



outside dim.	300 x 140 x 40 mm
inside dim.	280 x 130 x 27 mm



outside dim.	300 x 140 x 70 mm
inside dim.	280 x 130 x 55 mm



outside dim.	300 x 140 x 100 mm
inside dim.	280 x 130 x 84 mm

-  Y 410.07
-  Y 420.07
-  Y 430.07
-  Y 440.07
-  Y 450.07
-  Y 460.07

weight 567.5 gr

-  Y 411.07
-  Y 421.07
-  Y 431.07
-  Y 441.07
-  Y 451.07
-  Y 461.07

weight 630 gr

-  Y 410.10
-  Y 420.10
-  Y 430.10
-  Y 440.10
-  Y 450.10
-  Y 460.10


weight 625.5 gr

-  Y 411.10
-  Y 421.10
-  Y 431.10
-  Y 441.10
-  Y 451.10
-  Y 461.10

weight 625 gr



Bottom Non Perforated



Bottom Perforated



Lid Perforated



Y 400.01

Y 400.02

Y 400.03

Y 400.04

Y 400.05

Y 400.06

weight 325.5 gr

	outside dim.	inside dim.	weight
Y 405.04	300 x 140 x 40 mm	280 x 130 x 27 mm	190 gr
Y 405.07	300 x 140 x 70 mm	280 x 130 x 55 mm	241 gr
Y 405.10	300 x 140 x 100 mm	280 x 130 x 84 mm	305 gr

	outside dim.	inside dim.	weight
Y 406.04	300 x 140 x 40 mm	280 x 130 x 27 mm	258.5 gr
Y 406.07	300 x 140 x 70 mm	280 x 130 x 55 mm	304.5 gr
Y 406.10	300 x 140 x 100 mm	280 x 130 x 84 mm	326.5 gr



300 x 140 mm

Special Accessories

Silicone Mat

Drape Retainer

Drape



	outside dim.	weight
Y 850.00	270 x 125 mm	121 gr

Y 495.00	280 x 125 mm
----------	--------------

Y 095.01	60 x 40 cm
----------	------------

Lid and Bottom Perforated



Types and Dimensions

Y 510.04

weight 783 gr



outside dim.	310 x 190 x 40 mm
inside dim.	280 x 182 x 35 mm

Y 510.06

weight 862 gr



outside dim.	310 x 190 x 65 mm
inside dim.	280 x 182 x 55 mm

Y 510.13

weight 1022 gr




outside dim.	310 x 190 x 130 mm
inside dim.	280 x 182 x 124 mm

Bottom Perforated



Lid Perforated



	outside dim.	inside dim.	weight
Y 506.04	310 x 190 x 40 mm	280 x 182 x 35 mm	341,5 gr
Y 506.06	310 x 190 x 65 mm	280 x 182 x 55 mm	415 gr
Y 506.13	310 x 190 x 130 mm	280 x 182 x 124 mm	570 gr



	outside dim.	weight
Y 500.01	310 x 190 mm	417,5 gr
Y 500.02	310 x 190 mm	424 gr


Special Accessories

Drape Retainer



Y 595.00 310 x 190 mm

Drape



Y 095.01	60 x 40 cm
Y 095.02	65 x 65 cm
Y 095.03	85 x 85 cm

Implant containers

Types and Dimensions

Lid Perforated - Bottom Perforated (Non Filter)



- Y 612.07
- Y 622.07
- Y 632.07
- Y 642.07
- Y 652.07
- Y 662.07

weight 654 gr

outside dim.	500 x 155 x 75 mm
inside dim.	495 x 155 x 75 mm



- Y 610.06
- Y 620.06
- Y 630.06
- Y 640.06
- Y 650.06
- Y 660.06


weight 368 gr

outside dim.	310 x 138 x 65 mm
inside dim.	300 x 138 x 65 mm

Scopy containers

Types and Dimensions

Lid Perforated - Bottom Perforated (Non Filter)



Y 600.02

weight 302.5 gr



outside dim.	450 x 70 x 70 mm
inside dim.	450 x 70 x 60 mm

Y 600.01

weight 152 gr



outside dim.	250 x 60 x 50 mm
inside dim.	250 x 55 x 50 mm



Removal containers

Types and Dimensions



outside dim.	285 x 280 x 100 mm
inside dim.	262 x 260 x 82 mm

	Y 315.10
	Y 325.10
	Y 335.10
	Y 345.10
	Y 355.10
	Y 365.10
weight	1677,5 gr



outside dim.	285 x 280 x 150 mm
inside dim.	262 x 260 x 132 mm


	Y 315.15
	Y 325.15
	Y 335.15
	Y 345.15
	Y 355.15
	Y 365.15
weight	1929 gr



outside dim.	580 x 280 x 150 mm
inside dim.	548 x 258 x 130 mm

	Y 115.15
	Y 125.15
	Y 135.15
	Y 145.15
	Y 155.15
	Y 165.15
weight	2883,5 gr

Bottom Non Perforated



Y 105.15	580 x 280 x 150 mm	2108 gr
----------	--------------------	---------



Y 305.10	285 x 280 x 150 mm	1096 gr
Y 305.15	285 x 280 x 150 mm	1327 gr

Removal containers

Lids - Non Perforated

Y 302.01

Y 302.02

Y 302.03

Y 302.04

Y 302.05

Y 302.06

weight 458 gr



285 x 280 mm

bottoms, lids and accessories



580 x 280 mm

Lids - Non Perforated

Y 102.01

Y 102.02




Y 102.03

Y 102.04

Y 102.05

Y 102.06

weight 875 gr

 Y 880.10 Y 880.12 Y 882.00

Y 882.00

Brush container. For 20 brushes

240 x 275 x 70 mm

Stainless steel

Brushes are taken out by rotating the handle under the container.

Containers are supplied together with wall mounting fixtures.



Y 880.11



Y 880.10 - 14

Brush container. For 20 brushes

240 x 275 x 70 mm

Aluminium

Brushes are taken out by rotating the handle under the container.


Containers are supplied together with wall mounting fixtures.



Y 880.13



Y 880.14



Gen. accessories





Y 800.00
Pack of 100
For 1/1, 3/4 and
Half Size Containers




Y 800.01
Single textile
For 1/1, 3/4 and
Half Size Containers



Y 802.00
Pack of 100
For Mini Containers
Y 804.00
100 Pcs/Package
For Dental Containers



Y 802.01
Single textile
For Mini Containers
Y 804.01
Single
For Dental Containers




Y 810.00
Pack of 100
Paper label with indicator
For 1/1, 3/4 and Half Size Containers



Y 812.00
Pack of 100
Paper label with indicator
For Mini and Dental Containers




Y 809.00 For Paper Filter
Y 809.01 For Textile Filter



Y 801.00
Filter Dispenser



Y 819.00
Pack of 100
Plastic Security Seal



Y 829.00
Dispenser for identification labels

Silicone Mat



Y 850.00	270 x 125 mm
Y 851.00	245 x 245 mm
Y 852.00	220 x 230 mm
Y 852.10	380 x 230 mm
Y 852.20	440 x 230 mm
Y 852.30	520 x 230 mm

Neurosurgery	Y 820.01
Spine	Y 820.02
ENT	Y 820.03
General	Y 820.04
Neurosurgery	Y 820.05
Micro	Y 820.06

Coloured identification labels
Printed to your specifications,
max. 16 letters.



Y 849.20
24X28X25
Ø 7-12 mm



Y 849.21
24X28X29
Ø 12-16 mm



Y 849.22
24X28X31
Ø 16-26 mm



Y 849.23
24X28X40
Ø 26-36 mm



Y 849.24
24X28X44
Ø 36-45 mm




Y 841.00 230 mm
Y 841.01 470 mm
Sliding Rails with
Securing Screws



Y 849.00
Fixation Clamp
50X20X12 mm



Y 849.01
Fixation Clamp
65X20X20 mm



Y 842.00 230 mm
Y 842.01 470mm
Limiting Bar with
Securing Screws

Y 830.01
Holding Pin
For Partition Sheets



Y 830.00
Holding Pin
For Wire Baskets



Y 845.00
Spiral Holding



Y 847.00
Double
Spiral Holding

Partition Sheets
Stainless Steel

Y 831.05 50X20 mm
Y 831.13 130X20 mm
Y 831.22 225X20 mm
Y 831.46 460X20 mm



Y 840.00
Divider, to be
used with
Y 841.00



Trolleys

Carrying



Y 890.00
Trolley
1157 x 706 x 1405 mm



Y 891.00
Trolley
1157 x 706 x 1405 mm

Num. index



Kat. No Cat. No	Sayfa Page	Kat. No Cat. No	Sayfa Page	Kat. No Cat. No	Sayfa Page	Kat. No Cat. No	Sayfa Page
Y 095.01	39-40	Y 113.20	15	Y 142.20	15	Y 200.01	22
Y 095.02	40	Y 113.26	15	Y 142.26	15	Y 200.02	22
Y 095.03	40	Y 115.15	42	Y 143.10	15	Y 200.03	22
Y 095.04	29-37	Y 120.10	14	Y 143.13	15	Y 200.04	22
Y 095.05	29-37	Y 120.13	14	Y 143.15	15	Y 200.05	22
Y 095.06	29	Y 120.15	14	Y 143.20	15	Y 200.06	22
Y 095.07	17-29	Y 120.20	14	Y 143.26	15	Y 203.01	22
Y 095.08	17-23	Y 120.26	14	Y 145.15	42	Y 203.02	22
Y 095.09	17	Y 120.62	32	Y 150.10	14	Y 203.03	22
Y 095.10	17	Y 120.68	32	Y 150.13	14	Y 203.04	22
Y 100.01	16	Y 121.10	14	Y 150.15	14	Y 203.05	22
Y 100.02	16	Y 121.13	14	Y 150.20	14	Y 203.06	22
Y 100.03	16	Y 121.15	14	Y 150.26	14	Y 205.10	22
Y 100.04	16	Y 121.20	14	Y 150.62	32	Y 205.13	22
Y 100.05	16	Y 121.26	14	Y 150.68	32	Y 205.15	22
Y 100.06	16	Y 121.62	32	Y 151.10	14	Y 206.10	22
Y 100.51	33	Y 121.68	32	Y 151.13	14	Y 206.13	22
Y 100.52	33	Y 122.10	15	Y 151.15	14	Y 206.15	22
Y 100.53	33	Y 122.13	15	Y 151.20	14	Y 210.10	20
Y 100.54	33	Y 122.15	15	Y 151.26	14	Y 210.13	20
Y 100.55	33	Y 122.20	15	Y 151.62	32	Y 210.15	20
Y 100.56	33	Y 122.26	15	Y 151.68	32	Y 211.10	20
Y 102.01	43	Y 123.10	15	Y 152.10	15	Y 211.13	20
Y 102.02	43	Y 123.13	15	Y 152.13	15	Y 211.15	20
Y 102.03	43	Y 123.15	15	Y 152.15	15	Y 212.10	21
Y 102.04	43	Y 123.20	15	Y 152.20	15	Y 212.13	21
Y 102.05	43	Y 123.26	15	Y 152.26	15	Y 212.15	21
Y 102.06	43	Y 125.15	42	Y 153.10	15	Y 213.10	21
Y 103.01	16	Y 130.10	14	Y 153.13	15	Y 213.13	21
Y 103.02	16	Y 130.13	14	Y 153.15	15	Y 213.15	21
Y 103.03	16	Y 130.15	14	Y 153.20	15	Y 220.10	20
Y 103.04	16	Y 130.20	14	Y 153.26	15	Y 220.13	20
Y 103.05	16	Y 130.26	14	Y 155.15	42	Y 220.15	20
Y 103.06	16	Y 130.62	32	Y 160.10	14	Y 221.10	20
Y 105.10	16	Y 130.68	32	Y 160.13	14	Y 221.13	20
Y 105.13	16	Y 131.10	14	Y 160.15	14	Y 221.15	20
Y 105.15	16-42	Y 131.13	14	Y 160.20	14	Y 222.10	21
Y 105.20	16	Y 131.15	14	Y 160.26	14	Y 222.13	21
Y 105.26	16	Y 131.20	14	Y 160.62	32	Y 222.15	21
Y 105.62	33	Y 131.26	14	Y 160.68	32	Y 223.10	21
Y 105.68	33	Y 131.62	32	Y 161.10	14	Y 223.13	21
Y 106.10	16	Y 131.68	32	Y 161.13	14	Y 223.15	21
Y 106.13	16	Y 132.10	15	Y 161.15	14	Y 230.10	20
Y 106.15	16	Y 132.13	15	Y 161.20	14	Y 230.13	20
Y 106.20	16	Y 132.15	15	Y 161.26	14	Y 230.15	20
Y 106.26	16	Y 132.20	15	Y 161.62	32	Y 231.10	20
Y 106.62	33	Y 132.26	15	Y 161.68	32	Y 231.13	20
Y 106.68	33	Y 133.10	15	Y 162.10	15	Y 231.15	20
Y 110.10	14	Y 133.13	15	Y 162.13	15	Y 232.10	21
Y 110.13	14	Y 133.15	15	Y 162.15	15	Y 232.13	21
Y 110.15	14	Y 133.20	15	Y 162.20	15	Y 232.15	21
Y 110.20	14	Y 133.26	15	Y 162.26	15	Y 233.10	21
Y 110.26	14	Y 135.15	42	Y 163.10	15	Y 233.13	21
Y 110.62	32	Y 140.10	14	Y 163.13	15	Y 233.15	21
Y 110.68	32	Y 140.13	14	Y 163.15	15	Y 240.10	20
Y 111.10	14	Y 140.15	14	Y 163.20	15	Y 240.13	20
Y 111.13	14	Y 140.20	14	Y 163.26	15	Y 240.15	20
Y 111.15	14	Y 140.26	14	Y 165.15	42	Y 241.10	20
Y 111.20	14	Y 140.62	32	Y 190.03	17	Y 241.13	20
Y 111.26	14	Y 140.68	32	Y 190.05	17	Y 241.15	20
Y 111.62	32	Y 141.10	14	Y 190.07	17	Y 242.10	21
Y 111.68	32	Y 141.13	14	Y 190.10	17	Y 242.13	21
Y 112.10	15	Y 141.15	14	Y 190.55	32	Y 242.15	21
Y 112.13	15	Y 141.20	14	Y 191.03	17	Y 243.10	21
Y 112.15	15	Y 141.26	14	Y 191.05	17	Y 243.13	21
Y 112.20	15	Y 141.62	32	Y 191.07	17	Y 243.15	21
Y 112.26	15	Y 141.68	32	Y 191.10	17	Y 250.10	20
Y 113.10	15	Y 142.10	15	Y 191.60	32	Y 250.13	20
Y 113.13	15	Y 142.13	15	Y 193.00	17	Y 250.15	20
Y 113.15	15	Y 142.15	15	Y 195.00	17	Y 251.10	20

Kat.No Cat.No	Sayfa Page	Kat.No Cat.No	Sayfa Page	Kat.No Cat.No	Sayfa Page	Kat.No Cat.No	Sayfa Page
Y 251.13	20	Y 311.08	36	Y 340.10	26	Y 363.26	27
Y 251.15	20	Y 311.10	26	Y 340.13	26	Y 365.10	42
Y 252.10	21	Y 311.13	26	Y 340.15	26	Y 365.15	42
Y 252.13	21	Y 311.15	26	Y 340.20	26	Y 380.03	29-37
Y 252.15	21	Y 311.20	26	Y 340.26	26	Y 390.05	29-37
Y 253.10	21	Y 311.26	26	Y 341.08	36	Y 390.07	29
Y 253.13	21	Y 312.10	27	Y 341.10	26	Y 390.10	29
Y 253.15	21	Y 312.13	27	Y 341.13	26	Y 393.00	29-37
Y 260.10	20	Y 312.15	27	Y 341.15	26	Y 395.00	29-37
Y 260.13	20	Y 312.20	27	Y 341.20	26	Y 400.01	39
Y 260.15	20	Y 312.26	27	Y 341.26	26	Y 400.02	39
Y 261.10	20	Y 313.10	27	Y 342.10	27	Y 400.03	39
Y 261.13	20	Y 313.13	27	Y 342.13	27	Y 400.04	39
Y 261.15	20	Y 313.15	27	Y 342.15	27	Y 400.05	39
Y 262.10	21	Y 313.20	27	Y 342.20	27	Y 400.06	39
Y 262.13	21	Y 313.26	27	Y 342.26	27	Y 405.04	39
Y 262.15	21	Y 315.10	42	Y 343.10	27	Y 405.07	39
Y 263.10	21	Y 315.15	42	Y 343.13	27	Y 405.10	39
Y 263.13	21	Y 320.05	36	Y 343.15	27	Y 406.04	39
Y 263.15	21	Y 320.08	36	Y 343.20	27	Y 406.07	39
Y 290.03	23	Y 320.10	26	Y 343.26	27	Y 406.10	39
Y 290.05	23	Y 320.13	26	Y 345.10	42	Y 410.04	38
Y 290.07	23	Y 320.15	26	Y 345.15	42	Y 410.07	38
Y 290.10	23	Y 320.20	26	Y 350.05	36	Y 410.10	38
Y 293.03	23	Y 320.26	26	Y 350.08	36	Y 411.04	38
Y 295.00	23	Y 321.08	36	Y 350.10	26	Y 411.07	38
Y 300.01	28	Y 321.10	26	Y 350.13	26	Y 411.10	38
Y 300.02	28	Y 321.13	26	Y 350.15	26	Y 420.04	38
Y 300.03	28	Y 321.15	26	Y 350.20	26	Y 420.07	38
Y 300.04	28	Y 321.20	26	Y 350.26	26	Y 420.10	38
Y 300.05	28	Y 321.26	26	Y 351.08	36	Y 421.04	38
Y 300.06	28	Y 322.10	27	Y 351.10	26	Y 421.07	38
Y 302.01	43	Y 322.13	27	Y 351.13	26	Y 421.10	38
Y 302.02	43	Y 322.15	27	Y 351.15	26	Y 430.04	38
Y 302.03	43	Y 322.20	27	Y 351.20	26	Y 430.07	38
Y 302.04	43	Y 322.26	27	Y 351.26	26	Y 430.10	38
Y 302.05	43	Y 323.10	27	Y 352.10	27	Y 431.04	38
Y 302.06	43	Y 323.13	27	Y 352.13	27	Y 431.07	38
Y 303.01	28	Y 323.15	27	Y 352.15	27	Y 431.10	38
Y 303.02	28	Y 323.20	27	Y 352.20	27	Y 440.04	38
Y 303.03	28	Y 323.26	27	Y 352.26	27	Y 440.07	38
Y 303.04	28	Y 325.10	42	Y 353.10	27	Y 440.10	38
Y 303.05	28	Y 325.15	42	Y 353.13	27	Y 441.04	38
Y 303.06	28	Y 330.05	36	Y 353.15	27	Y 441.07	38
Y 304.01	37	Y 330.08	36	Y 353.20	27	Y 441.10	38
Y 304.02	37	Y 330.10	26	Y 353.26	27	Y 450.04	38
Y 304.03	37	Y 330.13	26	Y 355.10	42	Y 450.07	38
Y 304.04	37	Y 330.15	26	Y 355.15	42	Y 450.10	38
Y 304.05	37	Y 330.20	26	Y 360.05	36	Y 451.04	38
Y 304.06	37	Y 330.26	26	Y 360.08	36	Y 451.07	38
Y 305.05	37	Y 331.08	36	Y 360.10	26	Y 451.10	38
Y 305.08	37	Y 331.10	26	Y 360.13	26	Y 460.04	38
Y 305.10	28-42	Y 331.13	26	Y 360.15	26	Y 460.07	38
Y 305.13	28	Y 331.15	26	Y 360.20	26	Y 460.10	38
Y 305.15	28-42	Y 331.20	26	Y 360.26	26	Y 461.04	38
Y 305.20	28	Y 331.26	26	Y 361.08	36	Y 461.07	38
Y 305.26	28	Y 332.10	27	Y 361.10	26	Y 461.10	38
Y 306.08	37	Y 332.13	27	Y 361.13	26	Y 495.00	39
Y 306.10	28	Y 332.15	27	Y 361.15	26	Y 500.01	40
Y 306.13	28	Y 332.20	27	Y 361.20	26	Y 500.02	40
Y 306.15	28	Y 332.26	27	Y 361.26	26	Y 506.04	40
Y 306.20	28	Y 333.10	27	Y 362.10	27	Y 506.06	40
Y 306.26	28	Y 333.13	27	Y 362.13	27	Y 506.13	40
Y 310.05	36	Y 333.15	27	Y 362.15	27	Y 510.04	40
Y 310.08	36	Y 333.20	27	Y 362.20	27	Y 510.06	40
Y 310.10	26	Y 333.26	27	Y 362.26	27	Y 510.13	40
Y 310.13	26	Y 335.10	42	Y 363.10	27	Y 595.00	40
Y 310.15	26	Y 335.15	42	Y 363.13	27	Y 600.01	41
Y 310.20	26	Y 340.05	36	Y 363.15	27	Y 600.02	41
Y 310.26	26	Y 340.08	36	Y 363.20	27	Y 610.06	41

Kat. No Cat. No	Sayfa Page	Kat. No Cat. No	Sayfa Page	Kat. No Cat. No	Sayfa Page	Kat. No Cat. No	Sayfa Page
Y 612.07	41						
Y 620.06	41						
Y 622.07	41						
Y 630.06	41						
Y 632.07	41						
Y 640.06	41						
Y 642.07	41						
Y 650.06	41						
Y 652.07	41						
Y 660.06	41						
Y 662.07	41						
Y 800.00	48						
Y 800.01	48						
Y 801.00	48						
Y 802.00	48						
Y 802.01	48						
Y 804.00	48						
Y 804.01	48						
Y 809.00	48						
Y 809.01	48						
Y 810.00	48						
Y 812.00	48						
Y 819.00	48						
Y 820.01	48						
Y 820.02	48						
Y 820.03	48						
Y 820.04	48						
Y 820.05	48						
Y 820.06	48						
Y 829.00	48						
Y 830.00	49						
Y 830.01	49						
Y 831.05	49						
Y 831.13	49						
Y 831.22	49						
Y 831.46	49						
Y 840.00	49						
Y 841.00	49						
Y 841.01	49						
Y 842.00	49						
Y 842.01	49						
Y 845.00	49						
Y 847.00	49						
Y 849.00	49						
Y 849.01	49						
Y 849.20	49						
Y 849.21	49						
Y 849.22	49						
Y 849.23	49						
Y 849.24	49						
Y 850.00	48-39						
Y 851.00	48-37						
Y 852.00	48						
Y 852.10	48						
Y 852.20	48						
Y 852.30	48						
Y 880.00	45						
Y 880.01	45						
Y 880.02	44						
Y 880.03	45						
Y 880.04	45						
Y 882.00	44						
Y 890.00	50						
Y 891.00	50						





TIBBİ ALETLER SANAYİ VE TİCARET A.Ş

www.bahadir.com.tr

e-mail : info@bahadir.com.tr

Factory and Head Office

Organize Sanayi Bölgesi 55300 Samsun / TURKEY
Phone : +90(362)266 70 00 Fax : +90(362)266 70 05

BAHADIR- Catalog nr. 2004-005